
**06855**  
**INCLINOMETRO D'OSUALDO**  
**D'OSUALDO INCLINOMETER**

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<b>1 GENERAL INFORMATION</b> .....	<b>5</b>
1.1 MANUFACTURER.....	5
1.2 LIMITATIONS OF USE AND SPECIAL WARNINGS.....	5
1.3 SIDE EFFECTS AND CONTRAINDICATIONS.....	5
1.4 USERS.....	5
1.5 UNPACKING.....	5
1.6 ENVIRONMENT OF USE.....	5
1.7 STORING.....	5
1.8 MAINTENANCE.....	5
1.9 CLEANING.....	5
1.10 DISINFECTION.....	5
1.11 ADDITIONAL INFORMATION.....	5
<b>2 GENERAL WARNINGS</b> .....	<b>5</b>
2.1 SYMBOLS ON THE MANUAL.....	6
2.2 SYMBOLS ON LABELS.....	6
2.3 SYMBOLS ON THE PACKAGING.....	6
<b>3 06855 – INCLINOMETER D’OSUALDO</b> .....	<b>7</b>
3.1 INTENDED USE AND DESCRIPTION.....	7
3.2 MATERIALS.....	7
3.3 PACKAGE CONTENT.....	7
3.4 SETUP AND USE.....	7
<b>4 06830 - ILIAC CRESTS ANALYZER</b> .....	<b>8</b>
4.1 INTENDED USE AND DESCRIPTION.....	8
4.2 MATERIALS.....	8
4.3 PACKAGE CONTENT.....	8
4.4 SETUP AND USE.....	8
<b>5 05005 – DELTA LEG N</b> .....	<b>9</b>
5.1 INTENDED USE AND DESCRIPTION.....	9
5.2 MATERIALS.....	9
5.3 PACKAGE CONTENT.....	9
5.4 ADDITIONAL INFORMATION.....	9
5.5 SETUP AND USE.....	10
<b>6 06730 - SET OF SHIMS G</b> .....	<b>10</b>
6.1 INTENDED USE AND DESCRIPTION.....	10
6.2 MATERIALS.....	10
6.3 PACKAGE CONTENT.....	10
6.4 SETUP AND USE.....	10
<b>7 05003 – ARCOMETER D’OSUALDO</b> .....	<b>11</b>
7.1 INTENDED USE AND DESCRIPTION.....	11
7.2 MATERIALS.....	11
7.3 PACKAGE CONTENT.....	11
7.4 ADDITIONAL INFORMATION.....	11
7.5 SETUP AND USE.....	12
7.6 PROCEDURE OF MEASUREMENT.....	12
7.7 ARCOMETER USE.....	13
<b>8 05004 - TRUNK ASYMMETRY INDICATOR</b> .....	<b>15</b>

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## 1 GENERAL INFORMATION

This manual contains information for practical, correct and safe use of the device and is intended for reading by specialised personnel and the user of the product. It is recommended that you read the entire manual carefully before using the product.

If you have any doubts or need clarification, please contact your dealer who will be able to advise you correctly.

The importance of reading and understanding the user manual is highlighted on the product by the following symbol:

Follow the instructions for use

### 1.1 Manufacturer

CHINESPORT S.P.A, Via Croazia, 2-33100 Udine Italy

+39 0432621621

<https://www.chinesport.com>

[chinesport@chinesport.it](mailto:chinesport@chinesport.it)

The company produces in accordance with the quality standard UNI EN ISO 13485:2016

### 1.2 LIMITATIONS OF USE AND SPECIAL WARNINGS

The devices must not be used:

- Without having carried out all the foreseen adjustments and dedicated adjustments
- When use is contraindicated or generates serious inconvenience to the patient
- With tampering and / or modifications on the original product
- In an improper configuration or when there is a suspicion that it may cause harm to the patient
- In a different way from the intended use

### 1.3 SIDE EFFECTS AND CONTRAINDICATIONS

There are no side effects. Always use the equipment in a suitable way. Knowledge of this instruction booklet is recommended.

Before proceeding with use, it is necessary to evaluate the patient's psycho-physical condition.

### 1.4 USERS

Trained medical and technical staff

It is up to the specialist to judge the physical fitness of the patient for whom the product is intended to be used.

Use under operator supervision is always recommended.

### 1.5 UNPACKING

Remove the aids from the packaging without using the cutter to avoid damaging them and properly dispose of the packaging.

Any damage due to transport must be promptly reported to the carrier and to the supplier: damaged aids must not be used.

### 1.6 ENVIRONMENT OF USE

Clinical/hospital environment. Use on flat, dry, stable and obstacle-free surfaces. Not for outdoor use.

### 1.7 STORING

Room temperature	Humidity	Pressure
5°C ~ 60°C	1% ~80%	70kPa~106kPa

### 1.8 MAINTENANCE

Keep constantly clean

FREQUENCY	OPERATION
12 months	General cleaning Visual and functional inspection of the accessory

### 1.9 CLEANING

Gently dust and in case of slight stains, pass the surface with a damp cloth or with a solution of neutral detergent, rinsing with a little water.

Dry avoiding direct exposure to sunlight.

**NEVER USE AGGRESSIVE OR ALCOHOLIC PRODUCTS OR DETERGENTS!**  
These detergents damage plastic materials and screen prints, invalidating the product warranty

### 1.10 DISINFECTION

For surface disinfection we recommend the use of disinfectant products that do not contain alcohol, ammonia or other aggressive substances.

### 1.11 ADDITIONAL INFORMATION

Further information available on the website:

[www.chinesport.com](http://www.chinesport.com)

## 2 GENERAL WARNINGS

Always refer to this manual for proper use of the device.

- Store the device in an environment that complies with the labels on the packaging and the specifications in this manual.
- The manufacturer shall not be liable, to the fullest extent permitted by applicable law, for any direct or

indirect, special, incidental or consequential damages caused by:

- Wrong use of the device
- Improper use of the device and outside of its intended use
- Using the device in environments not covered by this manual
- Use with unsuitable patients
- Use without checking the status of the device as described in the relevant paragraph
- Incorrect maintenance or lack of maintenance
- Use with parts or accessories that are not compatible or not approved by the manufacturer
- Incorrect disposal or disposal is other than as described in this manual.

	Keep this side at the top
	Do not use sharp blades to open

## 2.1 SYMBOLS ON THE MANUAL

	Warning!
---	----------

## 2.2 SYMBOLS ON LABELS

	Follow the instructions for use
	CE certified
	Dispose of properly
	Manufacturer
	Production date
	Indoor use only


## 2.3 SYMBOLS ON THE PACKAGING

	Fragile
	Recyclable
	Keep dry

### 3 06855 – INCLINOMETER D'OSUALDO

#### 3.1 INTENDED USE AND DESCRIPTION

Instrument for the evaluation of the gibbus (rotation angle of the trunk in anteroflexion) in the patient suffering from scoliosis and the Cobb angle on radiographs.

It consists of a nearly rectangular Plexiglas element bearing a goniometric scale, at the center of which there is a stick which is free to rotate and bears a bubble; the stick's free end bears the reading point for the goniometric scale.

The longer side of the rectangle features a recess to facilitate application to the patient (in the case of the spinous processes being protruding).

#### 3.2 MATERIALS

Plastic, PMMA.

#### 3.3 PACKAGE CONTENT

Inclinometer, 1pc


#### 3.4 SETUP AND USE

The device does not require installation and is ready for use.



#### Angle of trunk rotation measurement

The Scoliosis Research Society defines the angle of trunk rotation as the angle between a horizontal plane and a plane crossing the most prominent point of the thoracic gibbus.



1. Position the patient as described in Adams' test or forward - bending test: bent trunk, extended knees, arms dangling forward
2. Position the instrument, gently laid on the patient's back, searching for the highest point of rotation (apex of the gibbus)
3. Level the rotating stick
4. Read the degree of rotation
5. Measure any other gibbus

N.B. The degree of inclination can be read after removing the instrument from the patient's back as well (holding the stick still between the thumb and forefinger), without the risk of losing the value.

#### **Cobb angle measuring technique by means of the inclinometer**

The instrument's edge is drawn near the upper rim sign (the most inclined vertebra), the stick is then rotated until it is level; the degree of inclination is read on the goniometric scale. Then measure the lower rim sign.

The Cobb angle is given by the sum of the two inclination angles of the rim signs.

## **4 06830 - ILIAC CRESTS ANALYZER**

### **4.1 INTENDED USE AND DESCRIPTION**

Manual instrument suitable for assessing the heterometry of the lower limbs in discharge, for supine patient.

### **4.2 MATERIALS**

Plastic, PMMA.

### **4.3 PACKAGE CONTENT**

Iliac Crest Analyzer, 1 pc

### **4.4 SETUP AND USE**

The device does not require installation and is ready for use.

Position the articulated arms on the patient's iliac crests so as to check their flatness using the level / bubble installed on the central body.

In the event that the instrument highlights the presence of misalignment of the iliac crests it is generally possible to apply elevations from the lower side in order to achieve the alignment of the iliac crests.

### **06730 - SET OF SHIMS G**

By means of this evaluation, it is possible to quantify the extent of the difference in height of the iliac crests. Attention, however, this data should not be assimilated to an assessment of the heterometry of the lower limbs, in fact a difference in level of the pelvis (and therefore of the iliac crests) can also derive from other factors concerning the spine or the same conformation of the pelvis or the lower limbs, both for skeletal and functional reasons.

What we collect with this tool is therefore a datum that finds meaning only in the context of a more complete examination of the locomotor system.

## 5 05005 – DELTA LEG N

### 5.1 INTENDED USE AND DESCRIPTION

Manual instrument suitable for measuring the heterometry of the lower limbs in discharge, for the supine patient.

### 5.2 MATERIALS

Plastic, PMMA.

### 5.3 PACKAGE CONTENT

Body delta leg N, 1pc  
Left footboard, 1 pc  
Footboard right foot, 1pc

### 5.4 ADDITIONAL INFORMATION

#### A few considerations on lower limb heterometry:

As the term indicates, by heterometry, we mean a difference in length of the lower limbs, whatever the cause; with the same meaning, we can find various terms: disparity, discrepancy, unequal length; they are all essentially synonyms.

However, it is a good idea to avoid the term "dysmetria", occasionally used by some authors, because it is used with a completely different meaning in the neurological field.

In Anglo-Saxon literature, heterometry is usually indicated by the expression "leg length discrepancy - LLD".

#### What does heterometry come from?

There are two distinct categories:

a) anatomical or structural heterometry - is an effective difference in the length of the skeletal segments (due to congenital or acquired causes);

b) functional heterometry - when the difference is tied to muscular contractures, axial deviations, subluxations, articular lassitude, asymmetry of the pelvis, etc.; factors that when present in the limb, determine an apparent shortening of it.

#### When is a heterometry important?

Scientific literature isn't in agreement on what value of heterometry is to be considered clinically relevant. For some authors, 5 mm are already significant while for others, 2 cm can be unimportant.

Generally, they are personal opinions based on anecdotal observations.

At the root of these uncertainties is first of all the indeterminateness of the clinical measure, tied to the ample error margin associated with evaluation methods that have been available up to now.

### Methods for evaluating heterometry that are available to date

#### 1.

X-ray examination: this represents the "gold-standard" and is essential for programming a surgical operation in the clinically relevant forms.

On the other hand, it presents obvious limitations of use tied to radiological risk and therefore is usually not used for the more modest forms, which are also the most common.

#### 2.

Methods of clinical observation (commonly used in specialized outpatient departments and in gymnasiums):

a) tape measure method: using a tape measure, different anatomical references can be used:

- SIAS - internal malleolus
- SIAS - external malleolus
- Navel - internal malleolus

With these measurements, the error (evaluated by repeating the measurement various times by the same observer or by two different observers) can exceed one centimeter and in some controlled studies, differences greater than 2 cm have been unacknowledged.

b) Evaluation in upright posture while observing the height of the:

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With these methods, the difference can be "estimated" or calculated by placing lifts until achieving compensation (checked with a level).

These methods too - although used at length - have ample uncertainty and only rarely have been subjected to checking the error range linked with every measuring procedure.



From these brief considerations, the fact emerges that with the common clinical evaluation methods, heterometry of 0,5 cm cannot be evaluated and even an apparent difference of 1 cm is to be cautiously evaluated before making even a simple clinical decision such as adopting a lift.

**To obviate the evaluation limits available, a product has been ideated: the DeltaLeg. It is a manual, non-invasive instrument suited to value draining lower limb heterometry.**


It is made of a bar on which two orthogonal platforms are applied: one is fixed and acts as a reference while the second is mobile along the longitudinal axis of the bar and is equipped with a pointer.

The pointer indicates the positive or negative numerical value of the heterometry on a millimetric scale located on the upper surface of the bar. The value "zero" corresponds to the reference platform

## 5.5 SETUP AND USE

Note: The instrument's stiff structure, the precise sliding system of the mobile platform and the presence of the millimetric scale allow taking reliable and rapid measurements of the differences in length of the lower limbs with an uncertainty margin of a few millimeters.

## 6 06730 - SET OF SHIMS G

### 6.1 INTENDED USE AND DESCRIPTION

Set of shims for postural correction. Measure 44..

### 6.2 MATERIALS

Plastic.

### 6.3 PACKAGE CONTENT

Shims 1 cm, 6pcs  
Shim 0.5cm, 1pc

### 6.4 SETUP AND USE

Installation is not necessary.

Place under the patient's foot to evaluate postural correction.

**1** Patient in supine position and well aligned (head aligned with nose pointing towards the zenith, torso aligned, lower limbs parallel with rotulas pointing towards the zenith).

**2** Ensure that the instrument's longitudinal axis is well - aligned with the legs.

**3** Bring the fixed platform close to the left foot and apply moderate pressure in order to ensure good adhesion of the entire sole and the heel in particular (NB. in the event of a retraction of the Achilles tendon, only functional heterometry can be measured).

**4** Bring the left mobile platform closer while exercising equal pressure.

**5** The difference on the millimetric scale can be read at this point or after having removed the instrument (paying attention not to change the ratios between the platforms).

### Recording the difference value

Considering that the right side is the reference point, the difference can be expressed with a single number; negative in the event the left side is shorter and positive in the event the left side is longer. For example: -2 heterometry means that the lower left limb is 2 cm shorter; +3 means that the left limb is 3 cm longer than the right.

## 7 05003 – ARCOMETER D'OSUALDO

### 7.1 INTENDED USE AND DESCRIPTION

Measurement of the curves of the spine on the sagittal plane (kyphosis and lordosis).

### 7.2 MATERIALS

PMMA, PE.

### 7.3 PACKAGE CONTENT

Arcometer, 1pc  
 Additioinal arm, 1pc  
 Cobb angles table, 1pc


### 7.4 ADDITIONAL INFORMATION

Dorsal kyphosis is a physiological curve, with a not-well defined normality range; moreover we do not have clear cut end points to identify its beginning and its end. (Boseker EH e al, J Pediatric Orthopaedics, 2000; Wenger DR, Frick S, Spine, 1999).

In the clinical practice an upper limit of normality for kyphos between 40° and 45° is usually accepted.

Lordosis normality values are even more vague, nor there is an agreement as what end vertebrae regards (e.g. lower limit of L5 or upper limit of S1). (VAL: agreement concerning the end vertebrae to consider).

A normal subject continuously modifies its kyphosis and lordosis as consequence of postural adjustments.

#### Kyphosis can be described by two values only:

- The angle comprised under the curve (Cobb angle).
- The radius of the circumference.

#### Current evaluation of kyphosis and lordosis relies on radiographs

On a radiograph taken in upright posture (usually arms forward), in a side projection, Cobb angle can be measured (end vertebrae D1-D12; L1-L5) with the same procedure usually adopted for scoliosis.

In order to avoid excessive exposure to X-rays, radiographs should be taken only if strictly necessary; as consequence they cannot be requested in every clinical evaluation or to verify the efficacy of a brace or the action of exercises.

#### Technique to build Cobb angle

(Shea K.G. et al 1998).

Error connected with the measurement of Cobb angle on radiographs.

It is usually accepted a reading error between 3°-5° (up to 10° in congenital scoliosis)

It depends on:

- the experience of the examiner;
- the quality of radiographs;
- the amount of the curve;
- the availability of technological devices.


**The measurement of Cobb angle on radiographs can be quicker and more precise by using an inclinometer.**

- Whittle MW, Evans M: Instrument for measuring the Cobb angle in scoliosis. The Lancet, 1979, Feb: 414).

- D'Osualdo F., Schierano S., Iannis M; The levelprotractor, a new simple instrument to measure back hump and Cobb angle; Eur Med Phys 2000; 36:191-196.



On the same curve, Cobb angle varies with the length of the arch, radius value does not.

Starting from the end of the XIX century several instruments, to measure kyphosis from the surface, have appeared in scientific literature; among them we can quote:

- 1 Neugebauer device;
- 2 Debrunner kyphometer;

3 Myrin Inclinator; 4 Spinal Pantograph; 5 Flexicurve

However, all these instruments have been little used in every day clinical practice for several reasons.

7.5 SETUP AND USE


The arcometer is a manual device made of a ruled bar bearing three perpendicular arms on it: the first one fixed at one end, the intermediate one is mobile on two axes and the third one mobile on one axis only.

The ends of the three arms identify three points through which a single

circumference can be drawn. Using the arcometer we can measure the chord and the rise. These two values allow us to calculate the radius of the curve and Cobb angle by using a two entry table.

**NB:** the arcometer works according to the trigonometric assioma that through three points only a single circumference can be drawn. For this purpose the segment of the spine is arbitrary approximated to a part (arch) of a circumference.

Arcometer schema with its working way: note that, in order to make calculations easier, the intermediate arm has to be placed into the median point between the lateral arms.

The validation of the instrument was made studying two aspects:

**Concordance X-rays** – arcometer: the mean difference between two series of measurements was 2.9° Cobb; the concordance between the two methods of measurement was good, without systematic errors.

The length of the chord can have an influence on the results since a shorter chord bears a risk of underestimating the angle.


**Concordance intra-observer:** the concordance between two series of measurements was good, without systematic errors: The mean of the differences resulted 0.1° Cobb (90% of the cases difference within ±10° Cobb; 82% of the cases within ±5° Cobb; the difference shows a small increase for kyphosis values between 44° and 63.5°).

Bibliography: D’Oswaldo F., Schierano S., Iannis M. Validation of clinical measurement of kyphosis with a simple instrument, the arcometer. Spine 1997; 22:408-13.

7.6 PROCEDURE OF MEASUREMENT

To measure the spine on the sagittal plane (Kyphosis and lordosis) the arcometer is placed with the external bars at the end of the segment: the length is the chord. The central bar is then fixed on the middle point of the length, then gently approached to the back: on the same bar one can read the value of the rise.


Cobb’s angle along with the radius value can be easily read on a two entry table.

$A = \arcsin H/2 \cdot R$   
 R = Radius  
 H = Chord  
 A = Cobb angle  
 F = Rise

**7.7 ARCOMETER USE**

**DEFINING DORSAL KYPHOSIS**

The upper and the lower limit of kyphosis are identified clinically and made visible through a dermographic pencil

The arcometer can be used to measure kyphosis (and lordosis) in different postures:

- Standing, spontaneous posture (standard posture)
- Standing, arms forward (comparison with X-rays)
- Self-correction while sitting
- Self-correction while prone
- While wearing the brace

NB: being the measurement easy, quick and without any adverse effect, it can be repeated as many times as needed.

The instrument comes with table to calculate the Cobb angle (current use to integrate the clinical examination).

		Rise											
		1	2	3	4	5	6	7	8	9	10	11	12
Rise	50	9	18	27	36	45	54	63	71	79	87	95	103
	60	8	15	23	30	38	45	53	60	67	74	81	87
Cobb angle	70	7	13	20	26	33	39	45	52	58	64	70	76
	80	6	11	17	23	29	34	40	45	51	56	62	67
Chord	90	5	10	15	20	25	30	35	40	45	50	55	60
	100	4	9	14	18	23	27	32	36	41	45	50	54
	110	4	8	12	17	21	25	29	33	37	41	45	49
	120	4	8	11	15	19	23	27	30	34	38	42	45
	130	4	7	11	14	18	21	25	28	32	35	38	42
	140	3	7	10	13	16	20	23	26	29	33	36	39
	150	3	6	9	12	15	18	21	24	27	30	33	36
	160	3	6	9	11	14	17	20	23	26	29	31	34
	170	3	5	8	11	13	16	19	22	24	27	29	32
	180	3	5	8	10	13	15	18	20	23	25	28	30
	190	2	5	7	10	12	14	17	19	22	24	26	29
	200	2	5	7	9	11	14	16	18	21	23	25	27
	210	2	4	7	9	11	13	15	17	20	22	24	26
	220	2	4	6	8	10	12	15	17	19	21	23	25
	230	2	4	6	8	10	12	14	16	18	20	22	24
	240	2	4	6	8	10	11	13	15	17	19	21	23
	250	2	4	5	7	9	11	13	15	16	18	20	22
	260	2	4	5	7	9	11	12	14	16	18	19	21
	270	2	3	5	7	8	10	12	14	15	17	19	20
	280	2	3	5	7	8	10	11	13	15	16	18	20
	290	2	3	5	6	8	9	11	13	14	16	17	19

The values of the Cobb angle are easily obtained crossing the values of the rise and of the chord given by the instrument.

In order to use the measurement for comparisons through time it is wise to register not only Cobb value but also:

- the chord length
  - the end vertebrae of the chosen segment (the kyphotic curve may not correspond to the entire dorsal segment)
- In a curve considered pathological for its location (e.g. dorso-lumbar), Cobb angle is particularly important for comparisons through time (more than for its absolute value).

**KYPHOSIS MEASUREMENT IN STANDING SPONTANEOUS POSTURE**

**SELF-CORRECTION WHILE SITTING**

- Third: remove the brace and take a second measurement in standard posture without changing the length of the chord and using the skin mark as a reference for the intermediate foot.
- Forth: take a third measurement while making a self-correction.
- Since the segment measured in-brace usually does not correspond to the entire kyphosis, it is advisable to use the radius table (to prevent any confusion with Cobb values).

The more the measurement taken in-brace approaches to self correction measurement, the better the brace is working.


### SELF-CORRECTION WHILE PRONE

### MEASUREMENT WHILE WEARING THE BRACE

In order to measure kyphosis within the brace one needs:

- to adapt the brace by making a 2 cm wide round hole in the dorsal pad;
- to adapt the instrument by changing the intermediate arm: the arm used for this purpose is provided with a smaller foot in order to enter the hole of the dorsal pad.

*Pad hole*

*Arcometer with intermediate arm provided with a smaller foot to enter the pad hole*

To measure the correction within the brace it is advisable to take 3 measurements as it follows:

- First measurement while wearing the brace (fixing the length of the chord pay attention not to invade lordosis).
- Second: mark with a dermographic pencil the point corresponding to the hole (or intermediate foot).

## 8 05004 - TRUNK ASYMMETRY INDICATOR

### 8.1 INTENDED USE AND DESCRIPTION

Instrument composed of an orthopedic technician caliber, with a plumb line applied on a sliding guide. The combination of these two elements makes it possible to describe the trunk asymmetries in quantitative terms.

### 8.2 MATERIALS

PMMA.

### 8.3 PACKAGE CONTENT

Trunk indicator, 1pc  
Plombe line, 1pc

### 8.4 SETUP AND USE

To measure the offset of a column point from a midline - for example the line interglutea - the following procedure must be followed:

- .1 Position the caliper on the patient's body at the height of the point of the column to be evaluated
- .2 Adjust the plumb line so that it passes through the medial reference (intergluteal line or support base)
- .3 Read the value of the size of the caliber and calculate that of the shortest segment
- .4 Based on the measurements just taken, use the attached table which will provide the percentage of translation of the anatomical part of interest

## 9 TORSION METER

### 9.1 INTENDED USE AND DESCRIPTION

D'Osualdo-Corazza's torsion meter is a manual tool for the assessment of the degree of rotation of the trunk in standing position.

### 9.2 MATERIALS

Plastic, alluminium

### 9.3 PACKAGE CONTENT

Torsion meter, 1pc

### 9.4 SETUP AND USE

The assessment of a subject with suspected scoliosis is carried out in a stand-up position.

The base of the instrument is placed near the sacred; in this way zero reference plane is defined, with respect to which the alignment of the trunk on the transverse plane at the desired level will be read. When the mobile element is in turn approached to the back, it will slip in its housing by a number of degrees corresponding to the rotation of the trunk. An index on the housing accurately indicates the degrees of rotation.

The tool is suitable for a diversified use:

- in stand-up position, it accurately identifies a rotation of the trunk, one of the most sensitive clinical signs for identifying true scoliosis;
- we can compare the rotation in static with the rotation in anteroflexion (scoliometer, inclinometer), data not always superimposable and in turn correlated, but inaccurately, with the vertebral torsion measured on the radiograph;
- in orthostatism or even seated we can verify the possible change of rotation with a stretching exercise (information on defect reducibility);
- with a stabilized pelvis (preferably sitting) we can evaluate the excursion of the trunk in an active rotation movement (range and symmetry).

## 10 05001 – FOOT INDICATOR

### 10.1 INTENDED USE AND DESCRIPTION

Foot length tool.

### 10.2 MATERIALS

PMMA.

### 10.3 PACKAGE CONTENT

Foot indicator, 1pc

### 10.4 SETUP AND USE

This is a simple but accurate tool which offers some data which is generally overlooked, the foot length. Its importance is given by the fact that the culmination moment of growth of the foot is one or two years before the culmination moment of the height, thus during the developing stage considerable changes can be seen in the relation between foot and height. A rather high value of this relation, over 16% indicates a phase of fast growth of the foot at the first stages of puberty. A rather low value of this relation, less than 14,5% indicates a culmination moment of growth already completed. This value can be useful especially during a first visit when the clinician does not have enough information about the growth stages of the young patient.

## 11 06810 – PLUMB LINE

### 11.1 INTENDED USE AND DESCRIPTION

Tool for the evaluation of postural deviations and scoliotic attitudes.

### 11.2 MATERIALS


Metal.

### 11.3 PACKAGE CONTENT

Plumb line, 1pc

### 11.4 SETUP AND USE

Position the plumb line as in the image, proceed with the evaluation.

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Algo goniometer, 1pc

### 12.4 SETUP AND USE

The intensity of the pain is represented along a scale, in which one end indicates absent pain, the other a very strong pain.

In the Wong-Baker scale, more suitable for children over 3-4 years of age, the graduations are represented with smiley faces, the expression of which, by analogy, indicates the intensity of the painful experience (nothing, little, enough, strong, very strong).

In the numerical VAS scale, suitable for children over 8 years and adults, the painful experience graduations are indicated with numbers, from 0 to 10, where 0 indicates no pain and 10 the worst possible pain, with all intermediate graduations.

There are two measuring instruments associated with the algometer. Their composition is functional for the execution of angular and linear measurements: the assemblage point forms a mobile protractor for the measurement of the joint ranges, their contiguous position forms a graduated line, while a fixed protractor is printed on one of the two parts.

The complete evaluation can also be carried out with the support of Chinesport podoscopes:

### 13 02049 - HINDFOOT PROTRACTOR

#### 13.1 INTENDED USE AND DESCRIPTION

Tool for assessing the alignment of the hindfoot.

#### 13.2 MATERIALS

PMMA.

#### 13.3 PACKAGE CONTENT

Goniometer, 1pc

#### 13.4 SETUP AND USE

The hand is oriented parallel to the hindfoot axis and the value of valgus or varus is immediately readable in degrees. The limits of normality cannot be indicated in an absolute way, since they vary in relation to age and loading methods (mono or bipodalic); however one can roughly indicate as a threshold of attention a valgus greater than 10 °, any varus value and significant asymmetries between the two sides. In doubtful cases, it is not so much the absolute value as the trend over time that can guide us for a physiological situation or not and for this purpose the instrument, due to its ease of use and precision, is particularly useful.

#### 02990 – PODOSCOPIO LUX 02992 - PODOCOLOR

#### 02991 - PODOLUX

## **14 WARRANTY**

### **14.1 GENERAL CONDITIONS**

All Chinesport products are warranted against defects in materials or workmanship for a period of 24 months from the date of sale of the product, except for any exclusions, limitations or conditions defined at the time of delivery of the product.

The warranty is not valid in case of improper use, tampering with the device, abuse or modification of the product or for any use or operation not explicitly mentioned in this manual.

The warranty is not valid if the device has not been correctly maintained and documented in accordance with this manual, or if the instructions regarding storage, cleaning and sanitation are not followed.

The manufacturer is not responsible for any damage or injury or any situation caused by incorrect installation or configuration of the device or using equipment that does not comply with the instructions in the installation, assembly and operating manuals.

The manufacturer does not guarantee its products against defects or damages in the presence of extraordinary conditions such as: natural disasters, unauthorised maintenance and repairs, improper power supply (where applicable), use of parts or components or accessories not original, shipping damage not directly managed by the manufacturer, lack of maintenance, obvious negligence on the part of the user or operator.

The warranty does not cover consumables, rechargeable batteries, and in general, all material subject to wear, failures caused by impacts, falls, misuse or improper, accidental events, damage caused by transport. If the equipment is tampered with, the warranty expires automatically.

### **14.2 REPAIRS UNDER WARRANTY**

In the case of a report of defects in materials or workmanship, the manufacturer assesses whether the defect is covered by warranty.

Warranty repairs must be expressly requested and are to be understood in our laboratory, subject to authorisation and with the issue of the return number.

For products sent in their original packaging, the return shipment will be made freight free.

For warranty repairs, a fiscal document is required where the date of purchase is within the warranty period (sales note, purchase invoice, fiscal receipt).

Labour costs for warranty repair (when the warranty conditions are valid) are borne by the manufacturer.

Repairing or replacing a product does not renew or extend the terms and expiration dates of the warranty.

### **14.3 OUT-OF-WARRANTY REPAIRS**

Non-warranty products can be repaired by the manufacturer by returning them after having been authorised by the technical assistance service. The costs of repair, including shipping, materials and labour, are to be understood as being borne by the customer or the retailer. The parts and components being repaired are considered to be covered by warranty for 24 months from the date of receipt of the repaired device.

### **14.4 NON-DEFECTIVE PRODUCTS**

In the event that the manufacturer does not find any malfunction or defect in the returned products, it is concluded that the product is not to be considered as defective. Shipping and device management costs will be charged to the customer or distributor.

### **14.5 HOME REPAIRS**

In case of repair at the customer's premises, a written request must be made indicating the complete details of the applicant, the type of machine and the fault.

The kilometric cost for the technician's transfer is to be agreed in relation to the customer's urgency.

In the event that the machine in question is under warranty, only the costs of the transfer will be charged.

The time is counted from the departure of the technician from our laboratory until his return, the time of return will be estimated on the basis of the time spent on the outward journey.

### **14.6 SPARE PARTS**

A detailed list of all spare parts can be obtained from the manufacturer.

Spare parts are sold following a formal request for an offer to the technical assistance service. Processing times are related to the availability of the parts. Returns for spare parts are not accepted.

The payment will be cash on delivery unless otherwise agreed.

## 15 BIGLIOGRAPHY

### 15.1 Delta leg

Woerman AL, Binder-Macleod A. (1984) Leg Length discrepancy assessment: accuracy and precision in five clinical methods of evaluation. *The Journal of Orthopaedic & Sports Physical Therapy* 1984, vol 5(5); 230-239.

Beattie P, Isaacson K, Riddle DL, Rothstein J. (1990) Validity of derived measurements of leg-length differences obtained by use of a tape measure. *Physical Therapy* 1990 vol 70(3); 150-157.

Brady RJ, Skinner TM. & Gross MT. (2003) Limb Length Inequality: Clinical implications for assessment and intervention. *Journal of Orthopaedic & Sport Physical Therapy* 2003; 33(5); 221-234.

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D'Oswaldo F., Schierano S., Jannis M. Validation of clinical measurement of kyphosis with a simple instrument named arcometer. *Spine*, vol 22, n° 4 (408-413), Feb 1997

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D'Oswaldo F., Schierano S., Soldano F. M., Isola M. A new tridimensional approach to the evaluation of the spine through surface measurement: the BACES System. *Medical Engineering and Technology* 2002, vol 26, n° 3 (95-105)

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D'Oswaldo F., Schierano S., Cisotti C., Battel M.A., Saccavini M. Short- Term Effect Of Standard Milwaukee And Thoracolumbar Braces For Dorsal Kyphosis Evaluated With A Computerized Surface Technique (Italian Acronym: B.A.C.E.S.). *ESMAC Seventh annual Meeting – Leuven 2002; Gait & Posture* 2002, vol 16, Suppl 1 pp 150

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