SERENA - LZ112
NEVIO - LZ212
LETTINI DA VISITA
VISIT TABLE

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LIST OF CONTENTS

1.	GENERAL INFORMATION	12	9.	DISPOSAL.....	16
1.1	Manufacturer	12	10.	TECHNICAL SPECIFICATIONS	17
2.	PRODUCT DESCRIPTION AND INFORMATION	12	11.	TROUBLESHOOTING.....	17
2.1	DESCRIPTION	12	12.	SPARE PARTS.....	17
2.2	INTENDED USE	12	13.	WARRANTY.....	18
2.3	USERS.....	12	13.1	GENERAL CONDITIONS.....	18
2.4	ENVIRONMENT OF USE	12	13.2	REPAIRS UNDER WARRANTY	18
2.5	STORING	12	13.3	OUT-OF-WARRANTY REPAIRS	18
2.6	MANUFACTURER'S DECLARATIONS	12	13.4	NON-DEFECTIVE PRODUCTS.....	18
3.	GENERAL WARNINGS.....	12	13.5	HOME REPAIRS.....	18
3.1	SYMBOLS ON LABELS	13	13.6	SPARE PARTS	18
3.2	SYMBOLS IN THE MANUAL.....	13	14.	RECORDS OF OPERATIONS AND MAINTENANCE.....	19
3.3	SYMBOLS ON THE PACKAGING	13			
4.	THE CONTENT OF THE PACKAGING AND THE CHARACTERISTICS OF THE PRODUCT	14			
4.1	CONTENTS OF THE PACKAGE	14			
4.2	PRODUCT.....	14			
4.2.1	Model LZ1	14			
4.2.2	Model LZ2.....	14			
4.2.3	Power supply	14			
4.2.4	Hand control.....	14			
5.	ASSEMBLY AND FIRST USE	14			
5.1	ASSEMBLY.....	14			
5.2	FIRST USE.....	14			
6.	USE	15			
6.1	WARNINGS BEFORE USE	15			
6.2	USE OF THE DEVICE.....	15			
6.2.1	HEIGHT ADJUSTMENT.....	15			
6.2.2	HEADER ADJUSTMENT.....	15			
6.2.3	MOVEMENT OF THE EXAMINATION COUCH	15			
7.	MAINTENANCE AND SERVICE LIFE	16			
7.1	ROUTINE MAINTENANCE	16			
7.2	MAJFUNCTIONS, EXTRAORDINARY MAINTENANCE AND REPAIRS.....	16			
7.3	USEFUL LIFE	16			
8.	CLEANING	16			
8.1	WASHING.....	16			
8.2	DISINFECTION	16			

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PROSSIMA
VERIFICA

1. GENERAL INFORMATION

This manual contains information for practical, correct and safe use of the device and is intended for reading by specialised personnel and the user of the product. It is recommended that you read the entire manual carefully before using the product.

If you have any doubts or need clarification, please contact your dealer who will be able to advise you correctly.

The importance of reading and understanding the user manual is highlighted on the product by the following symbol:

Follow the instructions for use

1.1 Manufacturer

CHINESPORT S.P.A, Via Croazia, 2-33100 Udine Italy

+39 0432621621

<https://www.chinesport.com>

chinesport@chinesport.it

The company produces in accordance with the quality standard UNI EN ISO 13485:2016

2. PRODUCT DESCRIPTION AND INFORMATION

2.1 DESCRIPTION

Two-section examination couch with electric height adjustment.

The head section can be tilted manually by means of a rack and pinion system.

The examination couch is equipped with two panels covered with imitation leather and wheel kit for easy movement.

2.2 INTENDED USE

Basic examination couch.

The use of the device for purposes other than those defined in this manual is prohibited. The manufacturer declines any responsibility for damage to persons or property resulting from improper use of the device or in any case other than that provided for in this manual.

The manufacturer reserves the right to make changes to the product and the manual without prior notice in order to improve its characteristics and performance.

2.3 USERS

Trained medical and technical staff

It is up to the specialist to judge the physical fitness of the patient for whom the product is intended to be used.

Use under operator supervision is always recommended.

2.4 ENVIRONMENT OF USE

Clinical/hospital environment. Use on flat, dry, stable and obstacle-free surfaces. Not for outdoor use.

2.5 STORING

Store in a dry place at room temperature. Avoid excessive pressure and contact with discolouring materials. Avoid excessive exposure to direct sunlight.

2.6 MANUFACTURER'S DECLARATIONS

The manufacturer declares that

- the device is not a measuring instrument
- the device is not intended for clinical investigation
- the device is not sterile and is not for single use
- for a correct functioning and for the safety of the user, it is necessary that the ordinary maintenance operations are carried out as described in the relevant paragraph.
- the device cannot be used for purposes other than those stated in this manual

3. GENERAL WARNINGS

Always refer to this manual for proper use of the device.

- Store the device in an environment that complies with the labels on the packaging and the specifications in this manual.
- The useful life of the product is 10 years in accordance with the correct execution of the ordinary maintenance operations provided for in this manual. It is strictly forbidden to use the device beyond its stated useful life. At the end of its useful life, it is possible to proceed as described in the relevant paragraph
- The manufacturer shall not be liable, to the fullest extent permitted by applicable law, for any direct or indirect, special, incidental or consequential damages caused by:
 - o Wrong use of the device
 - o Improper use of the device and outside of its intended use
 - o Use of the device beyond the stated service life in this manual
 - o Using the device in environments not covered by this manual
 - o Use with unsuitable patients
 - o The distraction of operators or incorrect application of commands and adjustments
 - o Use without checking the status of the device as described in the relevant paragraph
 - o Incorrect maintenance or lack of maintenance
 - o Use with parts or accessories that are not compatible or not approved by the manufacturer
 - o Incorrect disposal or disposal is other than as described in this manual.

The device is equipped with labels to draw attention to particular dangers such as:

Hand crushing danger

Foot crushing danger

Therefore, pay particular attention when carrying out operations in areas adjacent to these symbols.

3.1 SYMBOLS ON LABELS

	Follow the instructions for use
	Applied part type B
	CE certified
	Dispose of properly
	Manufacturer
	Production date
	Indoor use only
	Electrical safety class 2

3.2 SYMBOLS IN THE MANUAL

	Attention!
	Electrical risk
	Caution
	Note
	Suggestion

3.3 SYMBOLS ON THE PACKAGING

	Fragile
	Recyclable
	Keep dry
	Keep this side at the top
	Do not use sharp blades to open

4. THE CONTENT OF THE PACKAGING AND THE CHARACTERISTICS OF THE PRODUCT

The product is delivered in a suitable cardboard packaging so that it can be received intact and functional. To open the package and remove its contents, pay attention to the warnings and symbols on the package itself.

Dispose of the packaging and the resulting material in an appropriate manner and follow the instructions in the packaging and in this manual.


4.1 CONTENTS OF THE PACKAGE

The package contains

- Instruction manual
- Product already assembled
- Packaging material to be disposed of


4.2 PRODUCT

4.2.1 Model LZ1


- | | |
|---------------------------|--|
| 1 Hand control | 6 Wheels to move the examination couch |
| 2 Header section | 7 Power supply |
| 3 Foot section | 8 Linear actuator |
| 4 Sheet (optional) | 9 Wheel drive control |
| 5 Sheet holder (optional) | 10 Hand control support |

4.2.2 Model LZ2


- | | |
|---------------------------|--|
| 1 Hand control | 6 Wheels to move the examination couch |
| 2 Header section | 7 Power supply |
| 3 Foot section | 8 Linear actuator |
| 4 Sheet (optional) | 9 Wheel drive control |
| 5 Sheet holder (optional) | 10 Hand control support |

4.2.3 Power supply

- | |
|------------------------------|
| 1 Network presence indicator |
|------------------------------|

4.2.4 Hand control

- | |
|------------------------------|
| 1 Command confirmation light |
| 2 Ascent |
| 3 Descent |

5. ASSEMBLY AND FIRST USE

5.1 ASSEMBLY

The device is delivered assembled and ready to use.

5.2 FIRST USE

Before proceeding with the first use, it is necessary to read this instruction manual.

Make sure that the product corresponds to what is described in this manual and that the contents of the package coincide with what is described in the relevant paragraph.

Before using the product for the first time:

- Make sure there are no components of the device inside the packaging (screws, nuts, handwheels...)
- Make sure there are no obvious signs of damage or tampering
- Check the rotation and sliding of the wheels
- Make sure there are no obvious signs of loosening of the screws
- Connect the power supply unit to an electrical outlet and check the operation of the height adjustment control
- Use the height adjustment control to bring the device to its minimum height and then bring it back to its maximum height.

6. USE

6.1 WARNINGS BEFORE USE

Before each use make sure that:


The product does not show any obvious signs of tampering or damage
The product has been sanitised in the parts in contact with the user
The physical and clinical condition of the user has been assessed and found to be appropriate for the use of the device.
The operating environment is in accordance with the provisions of this instruction manual.
There are no particular hazards in the areas around the device (shelves, obstacles, flammable materials, etc. ...)
The device is correctly connected to the power supply

6.2 USE OF THE DEVICE

6.2.1 HEIGHT ADJUSTMENT

The device has a height adjustment system by means of a linear actuator controlled by a hand control. To adjust the examination couch, press and hold the up or down button until the desired height is reached.



- 1 Ascent button
- 2 Descent button
- 3 Indicator light

During ascent or descent, the indicator light remains on for as long as the control buttons are pressed.

6.2.2 HEADER ADJUSTMENT

Lift the header by acting on the short side of the examination couch until you reach the desired position


To return to the initial position, lift the header up to the maximum inclination and take it back to the low level.

6.2.3 MOVEMENT OF THE EXAMINATION COUCH

You can easily move the examination couch using the standard equipment of wheels mounted on the examination couch.

To activate the wheel kit, press the kit operating levers located on the head and foot sides of the examination couch.


Once the examination couch is in the desired position, deactivate the kit by operating the levers from the bottom to the top

Do not move the examination couch with the patient positioned on it.

Make sure the examination couch is disconnected from the power supply before moving it.

Make sure you have deactivated the wheel kit before the patient gets on the bed.

7. MAINTENANCE AND SERVICE LIFE

7.1 ROUTINE MAINTENANCE

Scheduled maintenance is required for proper operation of the device, proper use, patient and operator safety, and performance assurance of the device.

Ordinary maintenance work may be carried out by specialised operators or authorised technical personnel.

PERIODICITY	CATEGORY	VERIFICATION
12 months	Visual inspection of the device	<p>The integrity of upholstery panels</p> <p>Frame integrity and absence of desoldered elements</p> <p>Presence of the four wheels of the wheel kit drive levers</p> <p>Presence of the hand control and integrity</p> <p>Presence of bumpers on the base frame</p> <p>Presence of rack for manual header adjustment</p> <p>Absence of evident deformations on the frames</p> <p>No friction noise or abnormal vibration</p>
	Functional verification	<p>Checking that the height adjustment is working properly</p> <p>Verification of tightening</p> <p>Correct pivoting of the wheels</p> <p>Correct rolling of the wheels</p> <p>Operation of the wheel kit drive levers</p> <p>Operation of the rack for adjusting the header</p> <p>Executing a full load lifting cycle</p>

7.2 MALFUNCTIONS, EXTRAORDINARY MAINTENANCE AND REPAIRS

Extraordinary maintenance work may only be carried out by personnel authorised by the manufacturer. Otherwise any warranty conditions will be immediately terminated. The manufacturer is not responsible for any damage caused using this product.

Any malfunctions found by the user must be promptly reported to the distributor or directly to the manufacturer and inhibit the use of the device.

Repairs may only be carried out by technical personnel authorised by the manufacturer and may include the withdrawal of the device in order to carry out the necessary repairs.

7.3 USEFUL LIFE

The service life of the device is defined at the beginning of this manual.

At the end of its useful life, you can proceed in the following ways:

- 3) Dispose of the device as described in the paragraph "Disposal".

- 4) Require the manufacturer to recondition and recertify the device so that it can continue to be used

As stated in the paragraph "General warnings", the manufacturer declines all responsibility for the use of the device beyond the useful life established in this manual.

8. CLEANING

It is necessary to clean the device at the end of each use if the device is intended for different users.

8.1 WASHING

Remove dust from metal parts using mild detergents and drying immediately after washing.

Clean the rubber parts and the antibrachial supports with a neutral product with a damp cloth and immediately dry the surfaces.

8.2 DISINFECTION

For disinfection of the product use disinfection with a low chlorine content, such as AMUCHINA® 10% or equivalent solutions with a concentration of 0.1% sodium hypochlorite and the following spectrum of action:

- bactericide in the presence of interfering substances according to EN 1276:1997;
- a fungicide in the presence of interfering substances according to EN 1650:1997;
- bactericidal surface test according to EN 13697:2001;
- fungicide surface test according to EN 13697:2001;
- active on Salmonella Typhimurium according to EN 13697:2001;
- active on: HIV, HAV - HBV, HCV;
- virucide according to EN14476:2005;
- active on H1N1 influenza virus according to EN14476:2005.

For safe use, refer to the instruction leaflet included in the product.

9. DISPOSAL

The symbol on the label of the equipment indicates that the waste must be subject to "separate collection".

Therefore, the user must either hand over the waste to the separate waste collection centres set up by the local authorities or hand it over to the retailer against the purchase of new equipment of an equivalent type. Separate waste collection and subsequent treatment, recovery and disposal operations promote the production of equipment with recycled materials and limit any negative effects on the environment and health caused by improper waste management. Illegal disposal of the product by the user entails the application of the administrative sanctions referred to in Article 50 et seq. of Legislative Decree no. 22/1997.

10. TECHNICAL SPECIFICATIONS

	LZ1	LZ2
Number of sections	2	2
Safe working load (SWL)	200kg	200kg
Maximum Patient Weight	180kg	180kg
Empty weight	50kg	50kg
Total weight including SWL workload	250kg	250kg
Height adjustment	Elettrica	Elettrica
Minimum/maximum examination couch height	50-80cm	56-92 cm
Headrest	68x70 cm	68x70 cm
Footrest	124x70 cm	124x70 cm
Power supply	100-240V, 50-60Hz, Max 145VA	100-240V, 50-60Hz, Max 145VA
IP protection	IP66	IP66
Headrest regulation	Manuale	Manuale
Tilting headrest	45°	45°
Fire resistance	1 IM	1 IM
Padding	40mm	40mm
Padding density	30kg/m ³	30kg/m ³
Wheels	65mm	65mm
Use / duty cycle	10% On max 2min	10% On max 2min

11. TROUBLESHOOTING

PROBLEM	SOLUTION
The bed cannot be adjusted in height	<p>Make sure the power supply unit is connected to the mains and to the actuator. Make sure the hand control is connected to the actuator</p> <p>Make sure the indicator light on the power supply unit is on.</p> <p>Make sure the examination couch is not already at the end of the movement (all low or all high)</p> <p>Press the up or down button and check that the indicator light on the hand control is on.</p> <p>Make sure the load on the table does not exceed the maximum permitted values indicated in the manual</p> <p>If all of the above tests fail, please contact your service representative.</p>
Wheels don't roll	<p>Make sure you have correctly operated the wheel kit from both sides of the table</p>

Make sure there are no elements that hinder the rotation of the wheels

If all of the above tests fail, please contact your service representative.

Make sure there are no objects or dirt residues blocking the drive

Make sure the control is not already activated

Make sure the examination couch is not loaded or that there is no patient on the examination couch

If all of the above tests fail, please contact your service representative.

Make sure the examination couch is not loaded or that there is no patient on the examination couch

Make sure the rack is correctly fixed to the examination couch


If all of the above tests fail, please contact your service representative.

The wheel drive control is not activated

The header does not adjust

12. SPARE PARTS

Contact your dealer or the manufacturer's technical support service for information on spare parts. Only use original spare parts supplied by the manufacturer or by authorised distributors.

1	Header	7	Caps
2	Footer	8	Power supply
3	Caps	9	Hand control
4	Sheet holder	10	Actuator
5	Sheet	11	Rack and pinion
6	Wheels		

13. WARRANTY

13.1 GENERAL CONDITIONS

All Chinesport products are warranted against defects in materials or workmanship for a period of 24 months from the date of sale of the product, except for any exclusions, limitations or conditions defined at the time of delivery of the product.

The warranty is not valid in case of improper use, tampering with the device, abuse or modification of the product or for any use or operation not explicitly mentioned in this manual.

The warranty is not valid if the device has not been correctly maintained and documented in accordance with this manual, or if the instructions regarding storage, cleaning and sanitation are not followed.

The manufacturer is not responsible for any damage or injury or any situation caused by incorrect installation or configuration of the device or using equipment that does not comply with the instructions in the installation, assembly and operating manuals.

The manufacturer does not guarantee its products against defects or damages in the presence of extraordinary conditions such as: natural disasters, unauthorised maintenance and repairs, improper power supply (where applicable), use of parts or components or accessories not original, shipping damage not directly managed by the manufacturer, lack of maintenance, obvious negligence on the part of the user or operator.

The warranty does not cover consumables, rechargeable batteries, and in general, all material subject to wear, failures caused by impacts, falls, misuse or improper, accidental events, damage caused by transport. If the equipment is tampered with, the warranty expires automatically.

13.2 REPAIRS UNDER WARRANTY

In the case of a report of defects in materials or workmanship, the manufacturer assesses whether the defect is covered by warranty.

Warranty repairs must be expressly requested and are to be understood in our laboratory, subject to authorisation and with the issue of the return number.

For products sent in their original packaging, the return shipment will be made freight free.

For warranty repairs, a fiscal document is required where the date of purchase is within the warranty period (sales note, purchase invoice, fiscal receipt).

Labour costs for warranty repair (when the warranty conditions are valid) are borne by the manufacturer.

Repairing or replacing a product does not renew or extend the terms and expiration dates of the warranty.

13.3 OUT-OF-WARRANTY REPAIRS

Non-warranty products can be repaired by the manufacturer by returning them after having been authorised by the technical assistance service. The costs of repair, including shipping, materials and labour, are to be understood as being borne by the customer or the retailer. The parts and components being repaired are considered to be covered by warranty for 24 months from the date of receipt of the repaired device.

13.4 NON-DEFECTIVE PRODUCTS

In the event that the manufacturer does not find any malfunction or defect in the returned products, it is concluded that the product is not to be considered as defective. Shipping and device management costs will be charged to the customer or distributor.

13.5 HOME REPAIRS

In case of repair at the customer's premises, a written request must be made indicating the complete details of the applicant, the type of machine and the fault.

The kilometric cost for the technician's transfer is to be agreed in relation to the customer's urgency.

In the event that the machine in question is under warranty, only the costs of the transfer will be charged.

The time is counted from the departure of the technician from our laboratory until his return, the time of return will be estimated on the basis of the time spent on the outward journey.

13.6 SPARE PARTS

A detailed list of all spare parts can be obtained from the manufacturer.

Spare parts are sold following a formal request for an offer to the technical assistance service. Processing times are related to the availability of the parts. Returns for spare parts are not accepted.

The payment will be cash on delivery unless otherwise agreed.

14. RECORDS OF OPERATIONS AND MAINTENANCE

APPARATUS

DATE OF INSTALLATION

SERIAL NUMBER

DATE
INTERVENTION

OPERATIONS PERFORMED

TECHNICIAN
IN CHARGE

SIGNATURE

NEXT CHECK

MADE IN ITALY